



Coag **S**

WHOLE BLOOD INR POC INSTRUMENT

classical photo-optical principle
liquid phase coagulation reaction
high sensitive recombinant thromboplastin
precise sampling system
small, compact, portable design



Coag S

WHOLE BLOOD POINT-OF-CARE COAGULATION SYSTEM

Key features:

- Screening INR test for professional user
- Compact hand-size design
- Assay from one drop of capillary blood
- Outstanding correlation with INR results of classical PT
- Screening from capillary blood, venous blood or plasma
- Coag S POC instrument
- Single use safety lancet
- Specialized pipette for precise sampling and reagent dispensing
- Disposable reagent cartridge
- Measuring cuvette

Hardware:

- One tempered and covered measuring channel
- Liquid crystal colour touchscreen
- Computer communication ports for software upgrade and data transfer
- DC Power supply or battery operation
- Optional battery

Software:

- Innovative user interface
- QC feature
- More than 1000 data storage capacity
- Connectivity to Laboratory or Clinical information systems



EUROPE

Diagon Ltd - Hungary
1047 Budapest, Baross u. 48-52.
phone.: (+36 1) 369 6500
fax: (+36 1) 369 6301
e-mail: diagon@diagon.com
Diagon GmbH – Austria
GZO – Zukunftsweg 1,
Siegendorf, Burgenland 7011
phone: (+43 2687) 42090
e-mail: diagonmbh@diagon.com
www.diagon.com

SOUTH AMERICA

Diagon Ltda - Brasil
Av. Joaquim José Diniz, 787
CEP 31910-520
Bairro Fernão Dias
Belo Horizonte/MG
phone: (+55 31) 3488 6777
fax: (+55 31) 3488 6778
e-mail: diagon@diagon.com.br
www.diagon.com.br

ASIA

PT Inti Diagonama Selaras - Indonesia
Perkantoran Plaza Pasifik A3/59 -
Jl. Raya Boulevard Barat
Kelapa Gading Permal Jakarta 14250
phone.: (+62 21) 458 40649, 458 40650
fax: (+62 21) 458 45024
e-mail: dsu@centrin.net.id