

# GET IN TOUCH



## Phone

+36 30 260 1280



## Email

sales@diagon.com



## Address

1047 Budapest, Baross str.  
48-52, Hungary



## Website

www.diagon.com



## SINGLE TUBE MULTI-PLATFORM CONTROL PORTFOLIO

Our multi-parameter  
assayed controls for  
3-part and 5-part  
differential hematology  
instruments with  
reticulocyte and nRBC  
parameters



“

3 AND 5 DIFF + RET/NRBC

1 CONTROL FOR ALL

UP TO 30 DAYS OPEN

6 MONTHS CLOSED STABILITY

”



## CONTROL BLOOD

The future of CONTROL SYSTEMS

Sysmex | Abbott | Coulter | Horiba-ABX | Mindray | Siemens

and many more

# Diagon Control Blood Materials Cross Ref. Table

Version V2023-I	D-Check D	D-Check D Plus	D-Check SDP	D-Check RET-HR	D-Check SDP Plus	D-Cell 0	EQC 1+2	EQC Diff 1+2
<b>ABBOTT</b>								
Cell-Dyn 3500 / 3700								
<b>HORIBA-ABX</b>								
MICROS Series								
PENTRA 60 / 60C+ / ES60								
PENTRA 80 / XL80								
<b>SIEMENS</b>								
ADVIA 60								
ADVIA 560								
<b>DIAGON</b>								
D-Cell 30 / 60 / 60CS								
Cellagon 3								
D-Cell 5D								
CellaGON 5Pro								
<b>COULTER</b>								
Act Diff Series								
Act 5 Diff series								
Act 5 Diff AL								
<b>SYSMEX</b>								
KX 21 / 21N								
XP-300								
SX Series								
XT Series								
XE Series								
XN Series								
<b>MINDRAY</b>								
MINDRAY BC-20s / BC-21s / BC-30s / BC-31s <sup>a</sup>								
MINDRAY BC-1800 / 2800 / 2900 / 3000 Plus / 3200								
MINDRAY BC-5120 / 5130 / 5140 / 5150 / 5000								
MINDRAY BC-5100 / 5180 / 5300 / 5380								
MINDRAY BC-5000 Vet								
MINDRAY BC-5500								
MINDRAY BC-5600 / 5800								
MINDRAY BC-6000 / 6200								
MINDRAY BC-6600 / 6800 / 6800 / 6800Plus / 760R / 760S								
MINDRAY BC-2800 Vet								

Version V2023-I	D-Check D	D-Check D Plus	D-Check SDP	D-Check RET-HR	D-Check SDP Plus	D-Cell 0	EQC 1+2	EQC Diff 1+2
<b>Dymind</b>								
DH-31 / 33 / 36								
DH-36Vet								
DH-50 / 51 / 53 / 56								
DH-71 / 76								
<b>OTHER</b>								
ANALYTICON Hemolyzer 5								
DIATRON Abacus / Arcus / Minicell / Twincell								
DIATRON Abacus 5								
ERBA Elite 5								
Genrui KT-60 <sup>a</sup>								
HUMAN HumaCount 5L								
LINEAR Stel 3								
NIHON KOHDEN MEK-6318 / 6400 / 6410 / 6420								
NIHON KOHDEN MEK-7222								
RAYTO RT-7200 / RT-7600 / RT-7600S								
SFRI Hemix 5								
BIOTA VABIO-360								
BIOTA VABIO-560								
<sup>a</sup> Values determined using Diagon brand reagents; <sup>b</sup> Additional calibrators available on request; <sup>c</sup> available on request								
<b>Valid Assays Sheets</b>								
	V2020-D1	V2020-DP1	V2020-SD1	V2020-RET1	V2020-SDP1			
	Recommended for this instrument							
	Recommended, but for non-diff parameters							
	Available on request							

Updated: 28-01-2023

HAEMATOLOGY  
EXTENDED  
STABILITY  
CONTROLS

SINGLE TUBE MULTI-PLATFORM  
CONTROLS

READY FOR USE  
LIQUID CONTROLS

3 LEVELS OF  
CONCENTRATION

INDEPENDENT THIRD PARTY  
CONTROLS

